

Prevalon AirTAP® Patient Repositioning System

Competency skill checklist



Name:
Date:

Steps	Completed	Comments
1. Places AirTAP under patient by following hospital protocol and ensures the top is aligned to patient's shoulder, ensuring the patient is centered.		
2. Ensures that correct side of glide sheet is facing up (blue side up/blue with black strips down).		
Turning and positioning patient	Completed	Comments
3. Inserts Air Pump Nozzle Hose properly into Quick Connect Valve, and ensures it is secure		
4. Makes sure Air Pump Power Cord is safely routed under the bed, not posing a trip hazard and away from the bed wheels. Ensures the Air Pump Power Cord is not plugged into the outlet on the bed.		
5. Makes sure Air Pump Hose is safely routed alongside the bed rail.		
6. Powers on Air Pump to aide in Q2 turning, boost, and lateral transfers.		
7. Ensures proper number of caregivers are present to boost or laterally transfer patient.		
8.		
9. Properly uses side straps and wedges to offload sacrum by pushing wedges under the glide sheet.		
10. Places the wedges with the appropriate side facing up (black up, blue down).		
11. Properly places "Anchor Wedge with tail" under thighs of patient and uses second wedge for torso.		
12. Ensures there is a proper gap between wedges to ensure proper sacral offloading (hand check).		

Prevalon AirTAP® Patient Repositioning System

Competency skill checklist

13. Uses Boosting Straps or Handles to boost patient efficiently.		
Maintenance	Completed	Comments
14. Only uses Sage microclimate management body pad with the AirTAP system.		
15. Cleanses soiled glide sheet and wedges with warm, damp cloth if necessary.		
16. Replaces microclimate management body pad when soiled.		
17. Discards system when no longer required for patient or upon discharge.		

Employee signature

Observer signature